

DrugConfirm™ Advanced

Instant Urine Drug Test Cups
Training Guide

Watch the video on our Training and Certification page.



The information in this presentation is a general overview on using Confirm BioSciences DrugConfirm™ Advanced urine drug test cup; however, this can be applied to drug test dip cards and other drug test kits.

Product Training Contents:

- Drug test cup overview
- Collection procedure
- Testing procedure
- Reading results
- Sending presumptive positives for confirmation

DrugConfirm™ Advanced Technical Information

It is important to read the package insert before using this screening device.

SCREENING RESULTS:

- The DrugConfirm™ Advanced Cup is used for Screening Only. Positive results obtained with this device are presumptive.
- Additional testing is necessary to confirm the presumptive positive results. Positive results should be “confirmed” by an alternate method such as GC/MS (Gas Chromatography/Mass Spectrometry).
- Professional judgment should be applied to any drug of abuse test result, particularly with preliminary positive results.

*The DrugConfirm™ Advanced screening device is a rapid, one step screening test for the simultaneous qualitative detection of multiple drugs and drug metabolites in human urine. For in vitro diagnostic use only.

DrugConfirm™ Advanced Product Overview

FEATURES & BENEFITS

- Lot number and expiration date included on the test pouch
- 5-16 panel configurations available
- Easy-to-use clean-catch cup
- Built-in adulterants and temperature strip
- Up to 99% accuracy
- CLIA Waived, FDA 510(K) Cleared, DOT Approved, CE Marked options

DrugConfirm™ Advanced Cut Off Levels


DRUG NAME	ABBR	DETECTION PERIOD	CUTOFF
AMPHETAMINE	AMP	2-4 Days	1000 ng/mL
ALCOHOL	ALC	Up to 80 Hours	300 ng/mL /500 ng/mL
BARBITURATES	BAR	4-7 Days	300 ng/mL
BENZODIAZEPINES	BZO	3-7 Days	300 ng/mL
BUPRENORPHINE	BUP	1-3 Days	10 ng/mL
COCAINE	COC	2-4 Days	300 ng/mL
FENTANYL	FEN	1-3 Days	200 ng/mL
METHADONE METABOLITE	EDDP	3-5 Days	300 ng/mL
MARIJUANA	THC	1-30 Days	50 ng/mL
METHADONE	MTD	1-3 Days	50 ng/mL

DrugConfirm™ Advanced Cut Off Levels

DRUG NAME	ABBR	DETECTION PERIOD	CUTOFF
METHAMPHETAMINE	mAMP	3-5 Days	1000 ng/mL
MORPHINE	MOR	2-4 Days	300 ng/mL
ECSTASY/MOLLY	MDMA	1-3 Days	500 ng/mL
OPIATES	OPI	2-4 Days	2000 ng/mL
OXYCODONE	OXY	2-4 Days	100 ng/mL
PHENCYCLIDINE	PCP	7-14 Days	25 ng/mL
PROPOXYPHENE	PPX	1-2 Days	300 ng/mL
SYNTHETIC MARIJUANA	K2/K3	Up to 32 Hours	50 ng/mL
TRICYCLIC ANTIDEPRESSANTS	TCA	2-10 Days	1000 ng/mL
TRAMADOL	TRA	N/A	200 ng/mL

Collection Procedure:


1. Bring the test cup to room temperature (10°C to 30°C) prior to testing. The test should be performed by a certified laboratory technician.



2. After collecting urine specimen, secure the cap and place on a flat surface.

3. Wait 5 minutes.

4. Peel off front label and interpret colored lines. Do not read results after 5 minutes.



Reading Results:



Negative (-): A rose-pink band is visible in each Control region (C) and the Test region (T). The presence of a colored band indicates the concentration of the specific drug is below the detection limit.

Presumptive Positive (+): A rose-pink band is visible in each Control region (C). No color band appears in the specific Test region (T). The absence of a colored band indicates the concentration of the specific drug is above the minimum detection limit.

Invalid: If a color band is not visible in each of the Control regions (C) or a color band only appears in the Test region (T), the test is invalid. Discard the test and administer a new one.

*NOTE: Any visible line, even a faint line, indicates a negative result.

Shipping Specimens to the Lab


PLEASE CONTACT

**SS@CONFIRMBIOSCIENCES.COM FOR
INFORMATION ABOUT LAB CONFIRMATION
TESTING SERVICES THROUGH OUR PARTNER
LAB, CRL.**

SECURE COLLECTION CONTAINERS

URINE COLLECTION CUPS:

Double check to make sure the cup cap is tightened securely. Chain of Custody Forms (CCF) should be provided with your test cups.



ALL ON-SITE DEVICE SAMPLES & LAB URINE SAMPLES:

After labeling and sealing the specimen tightly, place it into the biohazard specimen bag with

For questions please contact:

info@confirmbiosciences.com

800-908-5603

***Note: Instruct the donor to initial the security seal. Collector should verify the info provided by the donor and validate that the specimen was collected correctly.**

